



# Avanti High Performance Centrifuges JXN Series & J Series

Versatility. Performance. Reliability.



**BRILLIANCE**  
*at every turn.*



# ADDRESS EVERY POINT OF YOUR WORKFLOW, OPTIMIZATION AND PREPARATION

## SPECIALTY LABWARE

Disposable HarvestLine System liners simplify the processing of large samples

- Convenient sample recovery and long term sample storage
- Enhances operator safety and efficiency
- Easy-to-use, low-cost sample containment

1 liter wide mouth polypropylene bottle features that maximize your workflow

- Achieve 9,000 rpm in the **JLA-9.1000** rotor
- AutoVent plug and polyphenylene sulfide cap provide greater chemical resistance



JLA-9.1000 rotors support 4 x 1 Liter polypropylene bottles fit with the HarvestLine System liners

*Eliminate steps, maximize time and minimize cross-contamination with less hands-on sample handling.*

## WORKFLOW OPTIMIZATION FROM BEGINNING TO END OF SAMPLE PREPARATION

### Cell Harvesting/ Pelleting




#### CELL HARVESTING

Maximize cell harvest volume. Large-volume processing rotor gives you efficiency and high throughput by processing up to 6 liters.



### Supernatant Retention




#### EXOSOME PREPARATION (COMBINATION)

Isolate microparticles with the gold standard of exosome preparation. Efficient and fast sample preparation solutions customized for each step in your specific exosome research workflow, including small volume rotors **JS-24.38** and **JA-25.50**.



### Cell Lysis

### Debris Clearing




#### CELL HARVESTING

Exceptional high-throughout performance for scaling up protein production workflow. Save time and increase efficiency of cell pelleting with the **JCF-Z** continuous flow rotor. Our rotor system processes large sample volumes capable of reaching up to 800 mL pellet size.



### Product Processing



#### VIRAL VECTOR PREPARATION

Achieve greater than 100,000 x g and purify your exosome by density gradient centrifugation using our **Avanti JXN 30** in combination with the **JS-24.15** or **JA-14.50\*** rotors.



\*Max speed: 35,050 x g.

*Brilliance at Every Turn*

# INCREASE PERFORMANCE AND RELIABILITY



## BROAD FUNCTIONALITY AND APPLICATIONS

- Beckman Counter rotor range adds versatility and complements ultra speed range
- We offer three rotors (JA-30.50Ti, JS-24.38, JS-24.15) that are capable of reaching speeds greater than 100,000 x g
- Easy loading with lower height and ergonomic rotor design

## SAMPLE PROTECTIONS WITH EXCEPTIONAL TEMPERATURE CONTROL

- Extensive profile of rotors that are made of a single piece of aluminum, which is better conductive of temperature
- Our JXN-30 centrifuge runs all rotors at the ideal 4° C
- Sample loading and unloading while rotor is integrated into the centrifuge system makes for a more efficient and less cumbersome and labor-intensive process

## PROTECT YOUR SAMPLE AND STANDARDIZE METHODS

- We have been manufacturing high performance centrifuges for nearly 50 years with a consistent rotor library to optimize backward compatibility

## MAXIMIZE LAB SAFETY

- Our centrifuges are equipped with pharmaceutical-grade sterilizing filters
- Rotor biocontainment protects lab personnel from labware failure



# FASTEST HIGH PERFORMANCE CENTRIFUGE IN THE MARKET



## Avanti JXN Series

### BUILT FOR VERSATILITY AND PERFORMANCE

To overcome the everyday challenges of a shared laboratory facility, choose one of the two Avanti high performance centrifuges.

- Widespread use in biologics production
- Compatible with an extensive rotor library, so you can choose from a wide variety of high performance rotors to cover applications from large-volume pelleting to high speed subcellular fractionation to cell separations
- Large touchscreen interface and additional multi-user support enhances control and usability

	Avanti JXN 30	Avanti JXN 26
<b>Working Capacity</b>	4.0 Liter	6.0 Liter
<b>Speed</b>	Up to 110,500 x g RCF	Up to 81,200 x g RCF

## Avanti J-26S XP

### SIMPLICITY WITH PERFORMANCE

Simple interface with access to a large rotor library with speeds and volumes to support applications for cells, proteins and harvesting.



*Brilliance at Every Turn*

# AVANTI JXN USER-FRIENDLY SOFTWARE



## INTUITIVE INTERFACE

- **Personalized user profile** with user login and admin control
- **Experiment Diagnostics** shows real-time run display of time, speed, and chamber temperature in graph
- **Visual color status** gives early warning of any protocol variations
- **Run History** records all run details, error messages, and it can be exported in Excel format

## REMOTE MONITORING

- Network the Avanti JXN to set up, monitor and control your centrifuge remotely
- Receive alerts and diagnostics remotely via email or SMS



## LARGE INTUITIVE TOUCHSCREEN

- Very clear and visible for users
- Makes setting and viewing parameters very easy
- Rotor code identification provides the option of setting rcf or rpm
- Independent overspeed control (windage & inertia check) ensures safe rotor operation
- 10 levels of acceleration and deceleration control



## MOBILE CONTROL

- **Mobilefuge mobile application** controls, schedules and monitors runs remotely

## BUILT-IN SAFETY AND EFFICIENCY WITH ROTOR INNOVATIONS








- **Dynamic Rotor Inertia Check (DRIC)** makes sure that rotors are not run beyond their maximum speed
- **Lightweight 6 x 1 L JLA-8.1000 rotor** offers the most efficient pelleting for large volumes, and easy handling and cleaning because of its modular design
- **Foot pedal for hands-free door opening** and an integrated surface work space optimize everyday use
- **Proven reliability** with robust design for consistent results and maximum uptime, year after year
- **Exceptional performance** meeting evolving application needs, from academic research to production facilities, with up to 110,500 x g with the JXN 30 and up to 6 liters capacity with the JXN 26
- **Versatile, high-throughput sample processing** from 6 x 1 L bottles, to 28 x 50 mL conical tubes, to 24 microplates per run, saving time and ensuring sample integrity with fast acceleration and deceleration rates, combined with smooth runs at all speeds

# ROTOR SPECIFICATIONS

Rotor	Capacity	Speed (RPM/RCF)	Application								
			P	SF	LH	CM	NA	C	V	B	
 <b>JA-30.50 Ti</b>	400 mL (8 x 50 mL)	30,000 108,860	✓	✓	✓	✓	✓	✓	✓	✓	✓
 <b>JA-25.50</b>	400 mL (8 x 50 mL)	25,000 75,600	✓	✓	✓	✓	✓	✓	✓	✓	✓
 <b>JA-25.15</b>	360 mL (24 x 15 mL)	25,000 74,200	✓	✓	✓			✓	✓	✓	✓
 <b>JA-21</b>	180 mL (18 x 10 mL)	21,000 50,400	✓	✓	✓			✓	✓	✓	✓
 <b>JA-20.1</b>	480 mL (32 x 15 mL)	20,000 51,500	✓	✓	✓			✓	✓	✓	✓
 <b>JA-20</b>	400 mL (8 x 50 mL)	20,000 48,400	✓	✓	✓			✓	✓	✓	✓
 <b>JA-18.1</b>	43.2 mL (24 x 18 mL)	18,000 42,100	✓	✓	✓			✓	✓	✓	✓
 <b>JA-18</b>	1,000 mL (10 x 100 mL)	18,000 47,900	✓	✓	✓			✓	✓	✓	✓
 <b>JA-17</b>	700 mL (14 x 50 mL)	17,000 39,800	✓	✓	✓			✓	✓	✓	✓
 <b>JLA-16.250</b>	1,500 mL (6 x 250 mL)	16,000 38,400	✓	✓	✓	✓	✓	✓	✓		✓
 <b>JA-14.50</b>	800 mL (16 x 50 mL)	14,000 35,000	✓	✓	✓	✓	✓	✓	✓	✓	✓
 <b>JA-14</b>	1,500 mL (6 x 250 mL)	14,000 30,100	✓	✓	✓	✓	✓	✓	✓	✓	✓
 <b>JLA-12.500</b>	3,000 mL (6 x 500 mL)	12,000 26,890	✓	✓	✓	✓	✓	✓	✓		✓
 <b>JA-12</b>	600 mL (12 x 50 mL)	12,000 23,200	✓		✓			✓	✓	✓	✓
 <b>JLA-10.500</b>	3,000 mL (6 x 500 mL)	10,000 18,500	✓	✓	✓	✓	✓	✓	✓		✓

Brilliance at Every Turn

Application

Rotor	Capacity	Speed (RPM/RCF)	P	SF	LH	CM	NA	C	V	B
 <b>JA-10</b>	3,000 mL (6 x 500 mL)	10,000 17,700	✓	✓	✓	✓	✓	✓		✓
 <b>JLA-9.1000</b>	4,000 mL (4 x 1000 mL)	9,000 16 800	✓		✓	✓	✓	✓	✓	✓
 <b>JLA-8.1000</b>	6,000 mL (6 x 1000 mL)	8,000 15 970	✓		✓	✓	✓	✓	✓	✓
 <b>JS-24.38</b>	231 mL (6 x 38.5 mL)	24,000 103 900	✓	✓	✓	✓	✓	✓	✓	✓
 <b>JS-13.1</b>	300 mL (6 x 50 mL)	13,000 26 500	✓	✓	✓		✓	✓	✓	✓
 <b>JS-5.3*</b>	2,000 mL (4 x 500 mL)	5,300 6130	✓	✓	✓	✓	✓	✓	✓	✓
 <b>JS-4.0*</b>	4,000 mL (4 x 1000 mL)	4,000 4,050	✓	✓	✓		✓	✓	✓	✓





**WORKFLOW**

**P** (Protein) **SF** (Subcellular Fractions) **LH** (Lysate/Tissue Homogenates)

**APPLICATIONS KEY**

**CM** (Clearing Media) **NA** (Nucleic Acid) **C** (Cells) **V** (Viruses) **B** (Blood)

\*are not compatible with JXN-30.

JCF-Z	Workflow application	Capacity	Avanti JXN-30	Avanti JXN-26
 <b>Std. Core</b>	Debris cleaning, large particle separation and cell pelleting	400 mL	20,000 39,900	
 <b>Large Core</b>		800 mL	20,000 39,900	
 <b>Small Core</b>		200 mL	20,000 39,900	
 <b>Zonal Core</b>		1,650 mL	20,000 39,900	

## COMPLIANCE AND TECHNICAL SUPPORT

- At Beckman Coulter Life Sciences, engineering, sales, support, training and service work together to offer comprehensive and extensive customer-focused products
- Expert service engineering team strives for “Fix It Right the First Time”
- Certifications of Compliance



Specifications	Avanti JXN-26	Avanti JXN-30	Avanti J-26S XP*
<b>Max. Speed/g-Force</b>	26,000 rpm/81,770 x g	30,000 rpm/110,500 x g	26,000 rpm/81,770 x g
<b>Max. Capacity</b>	6 Liters	4 Liters	6 Liters
<b>Speed Control</b>	1,000 to 10,000: ± 10 rpm 10,001 to max: ±0.1%	1,000 to 10,000:± 10 rpm 10,001 to max: ±0.1%	Low Speed ± 10 rpm of set speed High Speed ± 0.1% rpm of set speed
<b>Set Temperature</b>	-10° C to 40° C in 1° increments	-20° C to 40° C in 1° increments	-10°C to 40°C in 1° increments
<b>Temperature Accuracy</b>	±2° C of chamber temperature (after equilibration)	±2° C of chamber temperature (after equilibration)	±2° C of chamber temperature (after equilibration)
<b>Ambient Operating Range</b>	16° C to 38° C	16° C to 38° C	16°C to 38°C
<b>Accel/Decel Profiles</b>	11/12	11/12	2/3
<b>Heat Output</b>	5,120 BTU/hr (1.5 kW)	5,120 BTU/hr (1.5 kW)	6900 BTU/hr (2.0 kW)
<b>Drive Type/Cooling</b>	SR drive/Air-Cooled	SR drive/Air-Cooled	SR drive/Air-cooled
<b>Refrigeration</b>	Non-CFC, non-ozone depleting refrigerant	Non-CFC, non-ozone depleting refrigerant	Refrigerant R452A (HFO)
<b>Sound Level</b>	<62 dBA (1 m in front of the instrument, 1.5 m above the floor)	<62 dBA (1 m in front of the instrument, 1.5 m above the floor)	57 dBA (0.91 m / 3 ft from instrument at maximum speed)
<b>Dimensions</b>	71 W x 86 D x 86 H cm (28 W x 34 D x 34 H in)	71 W x 86 D x 86 H cm (28 W x 34 D x 34 H in)	71 W x 86 D x 86 H cm (28 W x 34 D x 34 H in)
<b>Weight</b>	290 kg (640 lb)	310 kg (680 lb)	290 kg (640 lbs)

Providing over 75 years of global leadership in centrifugation, Beckman Coulter Life Sciences designs, manufactures, sells, and services a complete line of centrifuge systems. By offering unique rotors and innovative bottles, tubes and accessories, coupled with advanced centrifugation software, Beckman Coulter delivers intelligent centrifugation solutions to laboratory science.



\*Non-IVD Centrifuge.

© 2023 Beckman Coulter, Inc. All rights reserved. Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For Beckman Coulter's worldwide office locations and phone numbers, please visit "Contact Us" at [beckman.com](https://www.beckman.com)

2023-GBL-EN-101732-v1