

OneTouch Tips

A truly universal tip that fits virtually any pipettor – OneTouch tips provide an innovative, ergonomic and environmentally friendly liquid handling solution



Innovative – Pipetting like you’ve never felt it
















A revolutionary tip that is manufactured from two resins, OneTouch tips provide a soft color-coded top collar and a rigid, ultra-clear polypropylene body.

Ergonomic – Soft tip requires less pressure to create an optimal seal

Tips reduce the necessary sealing and reciprocal ejection force felt with traditional single resin polypropylene tips. Tips are excellent for multichannel pipettors, and reduce pipettor wear.



OneTouch Tips

	Cat. No.	Vol. Range	Type	Graduations	Packaging	Pack Type	Quantity
10 µl Standard							
	10300	1-20 µl	Standard	1 µl, 5 µl, 10 µl	Non-Sterile, 96 Tips/Rack		10 Racks/Pack, 5 Packs/Case
	10310	1-20 µl	Standard	1 µl, 5 µl, 10 µl	Sterile, 96 Tips/Rack		10 Racks/Pack, 5 Packs/Case
10 µl Barrier							
	10320	1-20 µl	Barrier	1 µl, 5 µl, 10 µl	Sterile, 96 Tips/Rack		10 Racks/Pack, 5 Packs/Case
200 µl Standard							
	10330	1-330 µl	Standard	10 µl, 50 µl, 100 µl	Non-Sterile, 96 Tips/Rack		10 Racks/Pack, 5 Packs/Case
	10340	1-330 µl	Standard	10 µl, 50 µl, 100 µl	Sterile, 96 Tips/Rack		10 Racks/Pack, 5 Packs/Case
200 µl Barrier							
	10350	1-330 µl	Barrier	10 µl, 50 µl, 100 µl	Sterile, 96 Tips/Rack		10 Racks/Pack, 5 Packs/Case
1000 µl Standard							
	10360	50-1250 µl	Standard	200 µl, 500 µl, 1000 µl	Non-Sterile, 96 Tips/Rack		8 Racks/Pack, 5 Packs/Case
	10370	50-1250 µl	Standard	200 µl, 500 µl, 1000 µl	Sterile, 96 Tips/Rack		8 Racks/Pack, 5 Packs/Case
1000 µl Barrier							
	10380	50-1250 µl	Barrier	200 µl, 500 µl, 1000 µl	Sterile, 96 Tips/Rack		8 Racks/Pack, 5 Packs/Case